



3-D Rotating Mixer 18 rpm

LX910RM

3-D Rotating Mixer 18 rpm LX910RM

3-D Rotating Mixer 18 rpm LX910RM is a mixer with horizon, vertical and angled rotation with a speed of 18 rpm. Adjustable tube holder makes the product either horizontally or vertically or between the angles of horizon-vertical to rotate. Equipped with 2 adjustable rotating frames to meet different hybrid modes in one mixture.

Features

- ◆ Combination of adjustable rotating frame
- ◆ Variance of rotation modes - mild type to somersault
- ◆ Gentle press positive and negative rotation
- ◆ Provision with three types of rotating frame

Applications

Used in molecular hybridization, blood banking and various samples requiring mixing

Specifications

Model No	LX910RM
Speed range (fixed)	18 rpm
Voltage	AC 220 V / AC 110 V ; 50 / 60 Hz
Power	4 W
Dimension (W x D x H)	240 x 145 x 210 mm
Weight	1 kgs

Optional Accessories

Accessory No	Standard configuration	Description
RM01	1 pair (2 pieces)	5 centrifuge tubes 10 ml / 15 ml (\varnothing 15 mm - \varnothing 16 mm) + 6 centrifuge tubes 5 ml / 7 ml (\varnothing 12 mm - \varnothing 13 mm)
RM02	1 pair (2 pieces)	3 centrifuge tubes 50 ml + 18 centrifuge tubes 0.5 ml
RM03	1 pair (2 pieces)	18 centrifuge tubes 1.5 ml / 2.0 ml